

HIV1/2 Rapid Test

Wondfo One Step HIV1/2 Whole Blood/Serum/Plasma Test is a rapid immunochromatographic direct binding test for the qualitative detection of HIV1/2 antibodies in human whole blood, serum or plasma as an aid in the diagnosis of HIV infection.

▶ **Specimen**
Whole Blood/Serum/Plasma



▶ **Detection target**
HIV 1/2 antibodies

▶ **Storage temperature**
4~30°C

▶ **Shelf-life**
24 months



01 . Test Procedure



HIV infection is often diagnosed through rapid diagnostic tests (RDTs), which detect the presence or absence of HIV antibodies.

There were approximately 37.7 million people living with HIV at the end of 2020 with 1.7 million people becoming newly infected with HIV in 2021 globally. Sub-Saharan Africa is the most affected region, accounting for two-thirds of the global total of new HIV infections.

HIV infection is often diagnosed through rapid diagnostic tests (RDTs), which detect the presence or absence of HIV antibodies.

— Source: HIV data and statistics, WHO, 2020; <http://www.who.int/hiv/en/>

02 . Order Information

Item Code	Description	Pack size	Time to results
W50-C	One Step HIV 1/2 Whole Blood/Serum/Plasma Test	25T/kit	15 mins